2019-nCoV Antigentest (Lateral-Flow-Methode)

Katalog Nr.: **WG34P0016 (1T) / WG34P0017 (5T)**

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2019-nCoV Antigen Test (Lateral Flow Method)

Catalog No.: W634P0016 (1T) / W634P0017 (5T)

You must follow the test directions carefully to get an accurate result. The explanatory film is recommended to watch before using the test on: <https://haemato.de/schnelltest-antigen-selbsttest>

INTENDED USE

The Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an *in vitro* diagnostic test that is used for self-testing to detect the nucleocapsid protein of novel coronaviruses (2019-nCoV) antigen extracted from the nasal swab specimen. The test is intended as an aid in the diagnosis of coronavirus infection disease (COVID-19), which is caused by 2019-nCoV. For *in vitro* diagnostic use only. For self-testing use.

SUMMARY OF THE TEST

COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an immunochromatographic assay for rapid, qualitative detection of novel coronaviruses (2019-nCoV) antigen extracted from the nasal swab specimen. When specimen is dropped into cassette, capillary action carries the specimen to migrate along the membrane. When the 2019-nCoV antigen levels are equal to or above the target cutoff, there is a visible colored band in the "Test Region" (T), that indicates a positive result. Absence of this test colored band in the test region suggests a negative result. If the test is performed correctly, a line forms in the "Control Region" (C) of the test result window. The test provides preliminary test results. Negative results cannot exclude 2019-nCoV infection and they cannot be used as the sole basis for treatment or other management decision.

MAKE SURE YOUR TEST KIT CONTAINS:

- 2019-nCov Antigen Test Cassette
- Pre-filled Extraction Buffer
- Disposable Sterile Swab
- Biohazard Specimen Bag
- Procedure Card
- Instruction for Use

WHAT ELSE DO YOU NEED?

- Timer

WARNING AND PRECAUTION

1. Read the Instructions for use completely before using the product. Follow the instructions carefully. Failure to do so may result in inaccurate results.
2. Avoid getting the solutions contained within the kit in contact with your eyes, skin or clothing.
3. This kit is for external use, do not swallow.
4. Keep out of reach children.
5. The test kit is for single use, do not reuse any components of the test kit.
6. Do not use this test beyond the expiration date printed on the outer package. Always check expiration date prior to testing.
7. Do not touch the reaction area of the test cassette.
8. Do not use the kit if the sealed pouch is punctured or not well sealed.
9. DISPOSAL: All specimens and the used-kit bear an infectious risk. Discard all the test components in the provided biohazard specimen bag after use. The process of disposing the diagnostic kit must follow the local infectious disposal law.


STORAGE AND STABILITY

1. Store at 2~30°C in the sealed pouch up to the expiration date printed on the package. Do not freeze.
2. The test cassette should be used within 1 hour after taking out from the sealed pouch.
3. Keep away from sunlight, moisture and heat.
4. Kit contents are stable until the expiration date printed on the outer box.
5. The production date is printed on the outer box.


TEST PROCEDURE

Choose a location to do this test where you can run the test UNDISTURBED for 20 minutes. Wash and dry hands before you begin to perform the test. Bring the test components to room temperature (15~30°C). According to usability study of layman use which carried out by Wondfo. The test can be correctly performed for anyone age from 18 to 75. However, for the users age under 18, it is recommended to perform the test under adult supervision, while the users age over 75 should be assisted during the swabbing procedure.

1. Take out the Pre-filled Extraction Buffer Tube, unscrew the lid and place the tube in the tube rack (the tube rack is on or inside the test kit).
2. Take out the Test Cassette from the foil pouch and lay it out flat.



3. Remove the swab from the container, being careful not to touch the soft end, which is the absorbent tip.
4. Carefully insert the entire absorbent tip of the swab into a nostril (about 1.5~2 cm). Take samples from the nasal wall by firmly rotating the swab in a circular path five times against the nasal wall. Slowly remove the swab from the nostril. (This step should take approximately 15 seconds to ensure the collection of mucus and cells).
- Note:** Simply twirling the swab against one part of the inside of the nose or leaving the swab in the nose for 15 seconds is not a proper technique and may result in an insufficient sample.
5. Repeat the same process with the same swab in the other nostril.

***Caution: If the swab breaks during specimen collection, dispose the used testkit and repeat the steps with a new testkit.**

6. Insert the swab into the Pre-filled Extraction Buffer Tube and immerse the entire tip of swab into the Extraction Buffer. Rotate about 10 times.
7. Now hold the extraction tube with one hand and bend the handle of the swab down with the other hand until it breaks at the first break-point. Rescrew the lid and place the tube into the tube rack for 1 minute (use a timer).
8. Unscrew the small cap at the top of the Pre-filled Extraction Buffer Tube. Lay the Test Cassette flat and add 4 drops of processed specimen into the sample well.
9. Wait for 15~20 minutes and read the results. **Do not read results after 20 minutes.**

10. After the test is completed, place the Test Cassette in the Biohazard Specimen Bag and dispose all test kit materials via biohazard waste disposal protocol.
11. Re-apply hand sanitizer.

INTERPRETATION OF RESULT

Positive Result
Colored bands appear at both test line (T) and control line (C). This indicates a positive result for the 2019-nCoV antigen in the specimen. (See Q6)
It does not matter if one of the lines that make up the test line (T) is lighter or darker than other; the result is "Positive".

Negative Result
Colored band appears at control line (C) only. This indicates that the concentration of the 2019-nCoV antigen is zero or below the detection limit of the test. (See Q7)

Invalid Result
No visible colored band appears at control line after performing the test. The instructions may have not been followed correctly or the test may have deteriorated. It is recommended to repeat the sampling and testing procedures.

LIMITATIONS OF PROCEDURE

1. This reagent is designed to detect 2019-nCoV antigen in human nasal swab specimen.
2. Failure to follow the instructions for the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
4. This test has been developed for testing nasal swab specimen only.
5. The sample collection process will affect the accuracy of the test, such as improper sample collection, improper sample storage, or repeated freezing and thawing of the sample etc.
6. This reagent is a qualitative assay. As it is with any diagnostic procedure, a confirmed 2019-nCoV infection diagnosis should only be made by a physician after evaluating all clinical and laboratory findings.
7. Negative test results may occur if the level of antigen in a sample is below the detection limit of the test, or from improper sample collection, and the negative results are not intended to exclude other non 2019-nCoV virus infections.
8. Positive test results do not exclude co-infections with other pathogens and does not identify specific 2019-nCoV virus subtypes.
9. A negative test result does not rule out a coronavirus infection and does not exempt you from the applicable rules for spread control (e.g. contact restrictions and protective measures).

QUESTION & ANSWER

Q1. How does the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) work?

The Wondfo 2019-nCoV Antigen Test is a type of test called an antigen test. When you have COVID-19, the 2019-nCoV virus (the virus that causes COVID-19) can be present in your nasal secretions. The Wondfo 2019-nCoV Antigen Test can detect small parts of 2019-nCoV virus in your nasal secretions. These small parts of the 2019-nCoV virus are known as proteins or antigens.

Q2. What is the difference between a COVID-19 antigen, molecular, and antibody test?

There are different kinds of tests for diagnosing COVID-19. Molecular tests (also known as PCR tests) detect genetic material from the virus.

The Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) is an antigen test which detects small parts or proteins from the virus. Antigen tests are very specific for the virus but are not as sensitive as molecular tests.

Another type of test is an antibody test. A COVID-19 antibody test detects antibodies that have been made by your immune system in response to a previous COVID-19 infection.

Q3. Are the substances in this kit dangerous?

This kit does not contain measurable amounts of substances that could be dangerous to people. As a precaution, please carry out the following safety measures should you come into physical contact with the substances contained in the kit:

Inhalation:
If inhaled, move to fresh air. Should you develop a cough or difficulty breathing, seek medical attention from a doctor.

Eye contact:
Rinse the affected eye with plenty of tepid water for 15 min while separating the eyelids with the tips of your fingers. Remove contact lenses if you wear any. Should you experience redness or irritation of the affected eye, seek medical attention from a doctor.

Skin contact:
Wash the affected area with cool water and remove any contaminated clothing. Should you experience irritation of the skin or any other symptoms, seek medical attention from a doctor.

Ingestion:
Should you accidentally ingest components of the kit, gargle and avoid, if possible, vomiting. Contact the local poison control center or a doctor should you feel unsure or unwell.

Q4. Will this test hurt?

No, the disposal sterile swab is not sharp and it should not hurt. Sometimes the swab can feel slightly uncomfortable or tickly.

Q5. Why do I swab both nostrils?

Swabbing both nostrils gives you the best chance of collecting sufficient sample to generate an accurate result. It has been observed in some cases that only one nostril has detectable virus, so it is important to collect from both nostrils. Correct swabbing is important to obtain a correct result.

Q6. How accurate is the Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)?

In clinical study, 356 clinical case samples which include 138 confirmed as COVID-19 positive and 218 confirmed as COVID-19 negative by PCR assay were obtained for testing, and then compared the test results between Wondfo 2019-nCoV Antigen Test (Lateral Flow Method) and the PCR results. The results are shown below.

Reagent	PCR		Total
	Positive	Negative	
Wondfo 2019-nCoV Antigen Test (Lateral Flow Method)	135	2	137
	3	216	219
Total	138	218	356

Sensitivity: 97.83% (95%CI: 93.78%~99.55%)
Specificity: 99.08% (95%CI: 96.73%~99.89%)
Total agreement: 98.60% (95%CI: 96.75%~99.54%)

Q7. What does it mean if I have a positive test result?

A positive result means that you may have COVID-19 disease and it is important to be under the care of a healthcare provider, doctor or local clinical. It is likely you will be asked to isolate yourself at home to avoid spreading the virus to others, wear a face mask when recommended and wash your hands regularly with soap and water. A positive result does not in any way guarantee that you are or will be immune and therefore cannot (or can no longer) become infected. Do you also have cold symptoms in addition to the positive at-home test? Stay home for at least seven days and have as little contact as possible with others, including the people you live with. Use disposable tissues and throw them straight in the bin. Sneeze and cough into the crook of your elbow. Wash your hands regularly and wear a face mask. Are your symptoms getting worse (difficulty breathing, high fever, etc.)? Contact your doctor immediately.

Q8. What does it mean if I have a negative test result?

A negative result means the virus that causes COVID-19 was not found in your sample. A negative test result does not guarantee that you do not or have never had COVID-19, nor does it confirm whether or not you are currently contagious.

Do you have cold symptoms in addition to the negative at-home test? Since the at-home test does not provide complete certainty, you should assume that you have COVID-19. You can contact your doctor to find out if another test is needed. In the meantime, try to avoid leaving your home and have as little contact as possible with others, including the people you live with. Use disposable tissues and throw them straight in the bin. Sneeze and cough into the crook of your elbow. Wash your hands regularly and wear a face mask. Are your symptoms getting worse (difficulty breathing, high fever, etc.)? Contact your doctor/health provider immediately.

Q9. Can I get an incorrect "false" negative result with this test?

An incorrect or "false" negative can occur for any of the following reasons:

- Incorrectly reading test result as negative.
 - Not following the instructions for use carefully.
 - The amount of antigen in a sample may decrease as the duration of illness increases. As days post-symptom onset increase, antigen test results may be more likely to be negative compared to a molecular 2019-nCoV test. You may also test negative at the very beginning of your infection before you develop symptoms
- If this is the case, your healthcare provider, doctor or local clinical will consider the test result with all other aspects of your history such as symptoms and possible exposures to decide how to care for you. It is important you work with your healthcare provider to help you understand the next steps you should take.

Q10. Can I get an incorrect "false" positive result with this test?

An incorrect or "false" positive can occur for any of the following reasons:

- Incorrectly reading test result as positive.
 - Not following the instructions for use carefully.
- If this is the case, your healthcare provider, doctor or local clinical, will work with you to determine how best to care for you based on your test results along with medical history and your symptoms.

Q11. I have used the test but no colored band appears at control line (C). What should I do?

If there is no colored band appears at control line (C) within 15 minutes of performing the test, then the test has not worked. You should test again, using a new test, taking care to follow the instruction.

Q12. Can any medications or medical conditions affect the results?

The following substances show no interference with the test result of the test.

Type	Substance	Concentration
Allergic symptoms	Histamine Dihydrochloride	0.25 mg/mL
	Interferon alpha	1.25 ml/L
	Zanamivir	25 mg/L
	Ribavirin	375 mg/L
	Oseltamivir	187.5 mg/L
Antiviral drugs	Palamivir	750 mg/L
	Lopenavir	500 mg/L
	Ritonavir	125 mg/L
	Abidor	0.5 g/L
	Levofloxacin	1.25 g/L
	Azithromycin	2.5 g/L
Antibiotics	Ceftriaxone	2.5 g/L
	Meropenem	500 mg/L
	Tobramycin	0.25 ml/L

It is not known whether other drugs will affect the test result. If you have any questions, please contact your doctor or pharmacist and read the information leaflet that came with your medication before using the test.

Q13. What are the possible risks of this test?

- Possible Risks:
- Discomfort during the sampling
 - Incorrect test results (see Interpreting Results and Limitations Sections).

INDEX OF SYMBOL

IVD In Vitro Diagnostic Use	See Instruction for Use	Expiry Date	Tests per kit
Manufacturing Date	Keep Dry	LOT Batch Number	Store between 2~30°C
Keep away from Sunlight	Manufacturer	REF Catalog #	Do not reuse

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The BfArM Sonderzulassung only authorises the german Text of this Instruction for Use leaflet.

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HAEMATO PHARM GmbH